

To Whom It May Concern

Marketing Information Letter

Freka Connect ENFit/ProNeo, 1 ml (30VE)	9000788
Freka Connect ENFit/ProNeo 2,5 ml (30VE)	9000782
Freka Connect ENFit/ProNeo, 5 ml (30VE)	9000783
Freka Connect ENFit/ProNeo, 10 ml (10VE)	9000784XS
Freka Connect ENFit/ProNeo, 10 ml (30VE)	9000784S
Freka Connect ENFit/ProNeo,10 ml (100VE)	9000784
Freka Connect ENFit/ProNeo, 20 ml (10VE)	9000785XS
Freka Connect ENFit/ProNeo, 20 ml (30VE)	9000785S
Freka Connect ENFit/ProNeo,20 ml (100VE)	9000785
Freka Connect ENFit/ProNeo, 60 ml (10VE)	9000786XS
Freka Connect ENFit/ProNeo, 60 ml (30VE)	9000786S
Freka Connect ENFit/ProNeo, 60 ml (50VE)	9000786
Freka Connect ENFit/ProNeo, 100ml (10VE)	9000787XS
Freka Connect ENFit/ProNeo,100 ml (30VE)	9000787

Fresenius Kabi has launched Freka® Connect ENFit/ProNeo Enteral Syringes with the new ENFit connection acc. to ISO 80369-3. The Freka® Connect ENFit/ProNeo Enteral Syringes are available in different sizes and can be connected to:

- Transnasal and percutaneous Freka Feeding Tubes, ENFit
- Medication Port of enteral Giving Sets for pump and gravity assisted feeding
- The syringes are indicated for the administration of enteral nutrition and/or medication

The Freka® Connect ENFit/ProNeo Enteral Syringes are classified as single use – disposables Medical Devices. Single use can be stated as "One patient and one application".



**FRESENIUS
KABI**

Microbiological tests showed that the Freka® Connect ENFit/ProNeo Enteral Syringes can be used for up to 96 hours. For hygienic and functional reasons it is recommended to replace them at least every 96 hours.

The Freka® Connect ENFit/ProNeo Enteral Syringes have to be cleaned after each administration under running lukewarm water. Multiple flushing of the syringe body and the syringe piston is recommended. During feeding breaks they should be stored in a dry and clean environment and protected from direct sun light (please see also the label of the packaging).

For critical patients it is recommended to replace the Freka® Connect ENFit/ProNeo Enteral Syringes after each administration at least every 24 hours.

Bad Homburg, 03.05.2017

Sabine Krohn

Product Management
Medical Device Division
Fresenius Kabi AG

Heinrich Martens

Regulatory Affairs
Medical Device Division
Fresenius Kabi AG